Tiger Dental

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
 Achtung

Notizen/notes

A large grid of small dots arranged in approximately 25 rows and 40 columns, intended for taking notes.

English

Dear doctor, we are honored that you have chosen this innovative new product from Tiger Dental.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. For this reason, you may find further user notices online at tigerdental.com under the "Know-how" link. Our customer service is available for any additional questions.

The continual development of our products indicates our company's commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the "Downloads" link at tigerdental.com).

1. Manufacturer

Tiger Dental, Belruptstraße 59, 6900 Bregenz, Austria

2. Product description


The transpalatal arch is specially intended for use with the TopJet Distalizer. It is pre-shaped and transfers force to the molars on one or both sides. It is attached into the Goshgarian lock on the molar bands. U-shaped bends on both sides allow for a larger installation distance for the TopJet Distalizer.

3. Areas of application

Indications: The palatal arch takes up the force developed by the TopJet Distalizer. This force is passed to the two molars via the palatal arch. The palatal arch prevents rotation of the molars enabling frictionless development of force by the TopJet Distalizer. If required, expansion and/or derotation of the molars can be shaped into the palatal arch.

4. Details of the transpalatal arch

Size and description

The sections

- **A-segment:** This is where the TPA is glued into the Goshgarian lock. Retentions provide for a secure hold. Protrudes distally from the lock and can be pressed out again from a distal direction for easy removal of the TPA.
- **B-segment:** Responsible for the position of the TopJet clix to the centre of resistance of the molar.
- **C-segment:** Lengthens the distance to the mini anchor screw. The distance between the C-section to the mini-anchor screw should not be shorter than 14 mm.
- **D-segment:** This is where the T-Connector is fixed in place. Always adjust to the correct angle to the level of the occlusion! The distance between the D-section to the mini-anchor screw should not be shorter than 14 mm (measuring aid: angled micro-brush = 14 mm).
- **E-segment:** Returns the F-segment to the same line as the B-segment (sketch). Runs diagonally; the height and length can be adjusted here.
- **F-segment:** It is the only variable of the TPA and determines the total length. The length of the F-segment corresponds to the order size.
- **U-bend:** Formed by the C, D, and E segments.

5. Size selection

The interfossa distance is measured between the first two molars. The shape of the palate should also be taken into account in size selection. A high-arched/gothic palate may need the TPA to be a size bigger.


Size (F-Segment)	Interfossa distance
12 mm	42 mm
14 mm	44 mm
16 mm	46 mm
18 mm	48 mm
20 mm	50 mm
22 mm	52 mm
24 mm	54 mm

(e.g. interfossa distance 44 mm/F=14 mm: $15 + 14 + 15 = 44$ mm)


6. Bending instructions

6.1A Coffin lid bend: Creates the basic shape of the TPA. The F-segment should remain straight.

6.1B Expansion bend: The TPA should also expand the molars during distalization.

6.1A Coffin lid bend basic shape


6.1B Expansion bend: End by bending outwards

Action: Turn up segments A-E like a trapezium at the transition to segment F. The F-segment should remain straight.

Result: Ideal basic shape for TPA.

6.2 Comfort bend: This bend bridges a swelling in the palatal mucosa.


Comfort bend

Action: The B-segment is bent firmly outwards, so it is nearer to the horizontal.

Result: The side of the palatal arch avoids the mucosal swelling.

6.3 Position bend D-segment: This bend helps to bring the D-segment and thus the T-Connector into the right distance from the gingiva.


D-segment position bend possible in both directions

Action: Segments C and E are bent one after the other and to the same extent towards the middle.

Result: The D-segment and the T-Connector of the TopJet Distalizer gains more distance from the mucosa (if the reverse is done, less distance).

6.4 Anti-rotation bend A-segment: This bend in the A-segment sets up a counterrotation.

Action: Bend the A-segment to the middle.

Result: The medial rotation resulting from the distalization is compensated by this bend with lateral rotation of the molars.

7. Mounting the transpalatal arch

1. Bend the ends of the pre-formed TPA up like a trapezium at the transition E to F on both sides (coffin lid bend, Section 6.1A). This bend reflects the basic shape of the palate. The F-segment remains straight.
2. Bend segments B-C slightly outwards to gain more distance from the gingiva (comfort bend, Section 6.2).
3. Adjust segments C, D, E (approximately 1.5 mm distance from the gingiva) by turning the E-segment outwards or inwards. In each case, compensate for this bend by counter-turning the A-segment (positioning bend, Section 6.3).
5. If there are pressure points in the area of the T-Connector, bend the D-segment further away from the gums (position bend, Section 6.3).
6. The TPA should always be deployed gently for the desired distalization (expansion bend, Section 6.1B). Strongly rotated molars should if possible be rotated outwards before use of the TopJet Distalizer (Anti-rotation bend, Section 6.4).
7. At every patient check, the distance between the T-Connector and the gums (about 1.5 mm) must be measured with a probe. If necessary, alter the position of the bend.

8. Fastening in place

Use Triad® Gel to fasten the Topjet TPA in place; this is injected through the front and back opening of the Goshgarian lock (light curing). Whilst doing this, use a finger to press the TPA in a distal direction so that it is well positioned in the Goshgarian lock. Jiggle the assembly to test it.

9. Removal

Remove the Triad® Gel as best you can with the cutter and fissure drill. Then use the Weingart pliers to press the distal end of the TPA against the mesial end of the Goshgarian lock. Repeat on the other side. Loosen the TPA and remove on both sides.

Separating the TPA: If standard removal is not possible or the TPA is pressing deeply into the gum, it can be carefully separated using a cutting disc and removed in two parts. Insert the mini diamond disc into a blue anglepiece and with water cooling carefully separate the transpalatal arch. The TPA should always be secured with the Moskito clamp, needle holder or dental floss.


Emergency removal using separation

10. Notices for Single-use Products

The palatal arch is nearly ready to use out of the box, but is not sterile and should not be used more than once. The reprocessing of a once used palatal arch (recycling) and its reapplication to the patient is strictly prohibited.

11. Safety Instructions

This operating and installation manual contains tips and recommendations, but does not release you from personal responsibility! The use of the Palatal arch is allowed only by orthodontists and dentists according to RL 93/42/EWG. The user must comply with the provisions of the Medical Devices Directive and any other applicable laws and regulations.

 This product contains nickel.

12. Quality notice

Tiger Dental ensures the user a flawless product quality. The contents of this user manual are based on our own experience. The user is responsible for any manipulations made to the product. Tiger Dental shall not be liable for any manipulations made to their products by users in the absence of their influence.

Despite the fact that palatal arch a mature product, improvements and changes can be made that are not documented within the user manual. Additionally, we do not take responsibility for any printing errors.

13. Package contents

- Palate arch, 1 size, 5 pcs/pkg
- Palate arch assortment, 7 sizes (12–24 mm), 2 pcs per size 14 pcs/pkg

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