# Tiger Dental

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**Notizen / Notes**

A large grid of small dots for taking notes, consisting of 20 columns and 25 rows.

**English**

**Dear doctor,** we are honored that you have chosen this innovative new product from Tiger Dental.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. For this reason, you may find further user notices online at [tigerdental.com](http://tigerdental.com) under the “Know-how” link. Our customer service is available for any additional questions.

The continual development of our products indicates our company’s commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the “Downloads” link at [tigerdental.com](http://tigerdental.com)).

## 1. Manufacturer

Tiger Dental e.U., Belruptstrasse 59, 6900 Bregenz, Austria

## 2. Product Description

TopJet clix allows for reliable and incremental molar distalization through the use of a programmable telescoping arm. The operator benefits from the appliance’s optimized handling characteristics and patients benefit from its maximal comfort. The TopJet clix appliance offers appealing esthetics and enables compliance free one or two-sided distalization of the desired teeth of up to 0.8 mm per month. Due to its compact design, the installation of TopJet clix is amazingly simple and can be done in small and narrow palates. The TopJet has already been successfully used in patients 8 years and older. The appliance can be installed in just a few minutes. Add this feature to near invisibility and high patient comfort, and it is easy to see why this is such a great device. During the entire installation process, a thread secures the TopJet Distalizer to its packaging outside the mouth. This safety feature minimizes the risk of aspiration. The appliance can only be activated after the thread is cut. The appliance’s hinge-like connection to the palatal arch, called the T-connector, prevents distal tipping of the molars. Use of a transpalatal arch (TPA) ensures a lack of molar rotation and the TopJet’s friction-free power output. Operators can choose from either a 250 cN (black coding) or 360 cN (red coding) ver-

sion (see section “Selecting the TopJet clix”).

The TopJet clix can be special ordered for lengths in excess of 8 mm. At the maximum spring activation, the smallest distance between the micro implant and the TPA is 14 mm while the largest distance is 28 mm. This allows for a maximum distalization of 14 mm. The TopJet can be reactivated in 1.8 mm increments.

### Selection of TopJet clix

**Black Coding: TopJet clix 250 (250 cN)** Is mainly used in children to distalize the upper first molars only.

**Red Coding: TopJet clix 360 (360 cN)** For molars and premolar distalization. Can be deactivated to exert less force.

## 3. Field of Application

### Indications

1. Crowding in the maxillary anterior with sufficient space in the posterior for distalization (consider removal of third molars!)
2. Protruded maxillary anterior teeth
3. Lack of space in the maxillary buccal segments (i.e. peg laterals, missing/narrow primary teeth or impacted premolars).
4. Class II malocclusion with mesially migrated maxillary molars (also in skeletal class III/dental class II combination cases). The TopJet clix can be used without problems on one side or both sides.

### Contraindications

1. Periodontally compromised molars when using high force values
2. Short maxilla (distally directed second molars)

#### 4. Installation of TopJet clip

##### System Prerequisites


- 16/26 bands with windowed Goshgarian lingual sheaths
- TopJet transpalatal arch (prefabricated or using bending template)
- Dual Top micro implant (JetScrew) measuring: 2.0 mm x 10/12/14 mm, 2.5 mm x 12/14/16 mm and placed in the M4 position
- Off-angle handpiece with a hex driver attachment (113-MJ-204)
- A hooked explorer for activating and deactivating the TopJet
- Ligature cutter to cut the security thread
- TopJet-Pliers: to close the T- connector
- Tweezers: holding the TopJet clip during installation

##### Insertion in the mouth:


The TopJet must be checked for proper orientation, which is parallel to the occlusal plane prior to the placement. Depending on the obliquity of the axes of the screw, the orientation of the C-clip must be adjusted. To achieve this, the C-clip must be held flat with the Weingart-pliers and bent back or bent off at the bending point as needed, so that the T-Connector lies at the D-section of the transpalatal arch (TPA) tension-free. (Fig. 2-3 and 7-9). This bending ensures tension-free placement of the TopJet.

##### Clip-on: Pushing the C-clip onto the mini-anchor screw


Hold the TopJet with a fixation forceps and position the C-Clip on the screw head properly. **Clip-on:** Push the C-clip onto the head of the mini-anchor screw with gentle pressure using forceps (Weingart-pliers). **Caution:** The long end of the C-clip is orientated to median (Fig. 4-5). NB: The form of the C-clip mustn't be changed. Never bend the long end of the C-clip inwards and never change the diameter. Otherwise, the C-clip cannot be removed from the head of the mini-anchor screw.


1 The C-clip ensures easy clip-on / clip-off and secure hold of the TopJet Distalizer at the neck of the jet screw. The geometry of the C-clip must therefore not be altered!


2 Insertion angle of the screw

3 Bend the C-clip at the insertion angle

4 Clip-on the C-clip (before) at the Clip-on point


5 Clip-on (completet)


> Clip-off the C-Clips see Fig. 22

6 Insertion

7 Bending the C-Clip

8 Clip-on: Slide the clip (before)

9 Clip-on: C-clip fixed (completed)

10 Extend and position the T-connector

11 Close the T-connector

12 Coat the T-connector

13 Coat the screw head

## Connecting the T-connector

Now the T-connector is pulled out with a fixation forceps and positioned onto the TPA. The T-connector-pliers (modified tweed pliers) closes the brackets of the T-Connector so that it lies on the TPA with some space (Fig. 10-11).

With the mesial and distal connections engaged, the fixation or covering is made to increase comfort by applying light-curing resin. In this way, a rotational and angularly stable connection between the screw and the TopJet clix (Fig. 12-13) is created within seconds.

## Activation of TopJet clix

Now activate TopJet clix by cutting the security thread between the ball and the spring telescope. The patient will now feel a slight pressure now on the molars (Fig. 14-15).

## Easy Reactivation and Deactivation

Activation as well as deactivation are possible in seconds with TopJet clix using a hooked explorer or probe. During molar distalization, spring tension can be increased by hooking a probe behind the distal end of the telescope tube and pulling mesially. The spring will lock into place audibly. Too much activation? To deactivate the telescope simply disengage the spring lock with the same hooked explorer. When disengaging the spring lock, both tubes will slide distally, simultaneously decreasing the force on the molars. Fine adjustment in this fashion is extremely easy and can be selected at increments of 1.8 mm (Fig. 16-17).

14 Activation: Remove security thread

15 Insertion completed & activated

22 C-clip is now detached

23 Detaching the T-connector

16 Activation/reactivation

17 Deactivation

24 Swivel the Topjet upwards to loosen/detach the T-connector

25 Remove the Topjet

18 Remove Plastic-coat


19 Slide the c-clip off the screw


26 Traditional Removal: Press out the TPA from the distal side

27 In Case of Emergency: Cutting through the TPA

20 Clip-off the C-clip (before) at the Clip-off point

21 Clip-off (completet)

Topjet Distalizer full animation video:  
[www.topjet-distalizer.com](http://www.topjet-distalizer.com)

## 5. Removal of the TopJet clix distalizer & TPA

### Deactivation of the TopJet distalizer

First, the TopJet clix distalizer must be distally deactivated.

#### Clip-off: Removing the C-clip from the mini-anchor screw

Now remove the gel cover on the screw head with a cutter.

Clip-off: One part of the Weingart-pliers holds the long end of the C-clip at the clip-off point, the other one holds to the head of the mini-anchor screw. Now, push the C-clip against the mini-anchor screw with a little pressure (Fig. 18–22). It will also break the Triad® gel adhesive residues.

Caution: Do not remove the C-clip directly from the screw head by pulling the spring pin. Risk of loosening the screw!

#### Loosen the T-connector and remove the TopJet clix

After loosening the TopJet clix on the screw head, the spacer pin can be grasped with a Weingart-pliers. With a 45° swivel movement, the brackets of the T-connector open and it breaks out of the adhesive bond. The TopJet clix can be easily removed now (Fig. 23–25).

#### Removing a TPA in a classic way

In order to loosen the glued TPA from the Goshgarian lock, the slightly distally protruding part of the TPA is pressed out with a Weingart-pliers against the lock. (Fig. 26).

#### In Case of Emergency: Cutting Through the TPA

In the event that a traditional removal of the TPA is not possible (i.e. pressure from the TPA has ulcerated the palate), it is recommended to use the small diamond disc in an off-angle handpiece with water cooling to carefully sever the TPA. During this process ensure that the TPA is secured with mosquito pliers, needle-nosed pliers or with floss (Fig. 27).

## 6. Notices for Single-use Products

This product is nearly ready to use out of the box, but is not sterile and should not be used more than once. The reprocessing of this once used product (recycling) and its reapplication to the patient is strictly prohibited.

## 7. Safety Instructions

This operating and installation manual contains tips and recommendations, but does not release you from personal responsibility! The use of this product is allowed only by orthodontists and dentists according to MDD 93/42/EEC. The user must comply with the provisions of the Medical Devices Directive and any other applicable laws and regulations.

 This product contains nickel.


## 8. References of Quality

Tiger Dental assures the user a flawless quality of its products. The content of this manual is based on experience. The user is responsible for the processing of the products themselves. Tiger Dental is not liable for erroneous results or processing by the user that occur in the absence of influence of Tiger Dental. Although this product is a mature product, improvements and changes can be made, which may not be documented in the present manual. Similarly, we are not liable for printing errors.

## 9. Contents

REF 10664 TopJet clix 250  
REF 10665 TopJet clix 250+8  
REF 10666 TopJet clix 360  
REF 10667 TopJet clix 360+8


### 10. Technical specifications

Type	250	360	250+8	360+8
Max. strenght	250 cN	360 cN	250 cN	360 cN
Min. insertion lenght	14 mm	14 mm	22 mm	22 mm
Max. distal movement	14 mm	14 mm	14 mm	14 mm
Ratchet steps	5	5	5	5
Length per ratchet step	1.8 mm	1.8 mm	1.8 mm	1.8 mm

Ratchet step	mm	TJ 250	TJ 360
0	0	250 cN	360 cN
1	1,8	210 cN	300 cN
2	3,6	170 cN	240 cN
3	5,4	130 cN	180 cN
4	7,2	90 cN	120 cN
5	9,0	50 cN	60 cN

TopJet 250 = 40 grams per ratchet step, TopJet 360 = 60 grams per ratchet step

### 11. Description of the used symbols of ÖNORM EN ISO 15223-1

- Notified body
- Manufacturer
- Date of manufacture
- Lot / Batch number
- Reference / Article number
- Non sterile
- Do not reuse
- Consult instructions for use
- Caution

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